

Falco An alternative, more performant implementation of FastQC for high throughput sequence quality control(Galaxy Version 1.2.4+galaxy0)

Tool Parameters

Ran on the transcriptome

From webpage click on dataset and examine and compare the forward and reverse reads

fastp fast all-in-one preprocessing for FASTQ files(Galaxy Version 1.3.1+galaxy0)

Tool Parameters

Ran on the transcriptome pair-end collection

Looking at the html report looking at the number of reads quality versus quantity

HISAT2 A fast and sensitive alignment program(Galaxy Version 2.2.2+galaxy0)

Tool Parameters

Source for the reference genome

Use a genome from history

Built-in references were created using default options

Select the reference genome *

3: promicromonospora.fasta

s this a single or paired library

Paired-end Dataset Collection

Paired Collection *

146: fastp on collection 28: Read 1 output

Aligner

StringTie transcript assembly and quantification(Galaxy Version 2.2.3+galaxy0)

Tool Parameters

Input options

Short reads

Input short mapped reads *

219: HISAT2 on collection 192: aligned reads (BAM)

Aligned reads to transcript

StringTie merge transcripts(Galaxy Version 2.2.3+galaxy0)

Tool Parameters

Transcripts *

226: StringTie on collection 219: Assembled transcripts

switch to column select

GffCompare Compare assembled transcripts to a reference annotation (Galaxy Version 0.12.10+galaxy0)

Tool Parameters

GTF inputs for comparison *

233: StringTie merge on dataset 227-232 and collection 226

2: promicromonospora.gff

switch to column select

Compares to genome annotation

featureCounts Measure gene expression in RNA-Seq experiments from SAM or BAM files (Galaxy Version 2.1.1+galaxy0)

Tool Parameters

Alignment file *

219: HISAT2 on collection 192: aligned reads (BAM)

[The supplied input will be mapped over this tool.](#)

Specify strand information *

Unstranded

) (-s)

Gene annotation file

A GFF/GTF file in your history

Gene annotation file *

Does the input have read pairs?

Yes, paired-end and count them as 1 single fragment.